

Intensive Care Unit Solution



iCU complete

Intensive care unit complete

ELEGANZA 3XC
with Virtuoso mattress


iCU alive	Patient vital functions management	<ul style="list-style-type: none"> ✓ Smart CPR feature ✓ X-ray imaging 	
iCU safe	Prevention of pressure ulcers	<ul style="list-style-type: none"> ✓ Virtuoso alternating pressure mattress ✓ Ergoframe® positioning feature 	
	Prevention of respiratory complications	<ul style="list-style-type: none"> ✓ Cardiac Chair position 	
	Prevention of falls	<ul style="list-style-type: none"> ✓ Siderails with Protector ✓ Multi-zone bed exit alarm ✓ Brake alarm 	
iCU recovered	Accelerated recovery	<ul style="list-style-type: none"> ✓ Mobi-Lift® handle 	
iCU well	ICU tasks facilitation	<ul style="list-style-type: none"> ✓ Scales ✓ Foot controls with examination position 	

solution

ELEGANZA 5 with Virtuoso or OptiCare mattress

MULTICARE with Virtuoso or OptiCare mattress

- ✓ Smart CPR feature
- ✓ X-ray imaging



- ✓ Smart CPR feature
- ✓ X-ray imaging
- ✓ Automatic lateral therapy



- ✓ Virtuoso alternating pressure mattress
- ✓ OptiCare smart integrated mattress
- ✓ Ergoframe® positioning feature
- ✓ Lateral tilt



- ✓ Virtuoso alternating pressure mattress
- ✓ OptiCare smart integrated mattress
- ✓ Ergoframe® positioning feature
- ✓ Automatic lateral therapy



- ✓ Orthopnoeic Chair position
- ✓ 30°/45° backrest stop
- ✓ Lateral tilt 15°



- ✓ Orthopnoeic Chair position
- ✓ 30°/45° backrest stop
- ✓ Automatic lateral therapy



- ✓ Siderails with Protector
- ✓ Multi-zone bed exit alarm



- ✓ Siderails with Protector
- ✓ Multi-zone bed exit alarm
- ✓ i-Brake



- ✓ Mobi-Lift® handle + MobiRails + Lateral tilt for mobilisation



- ✓ Mobi-Lift® handle + MobiRails + Lateral tilt for mobilisation



- ✓ Scales
- ✓ Foot controls with examination position
- ✓ i-Drive Power
- ✓ Lateral tilt



- ✓ Scales
- ✓ Foot controls with examination position
- ✓ i-Drive Power
- ✓ Lateral tilt



icu **ALIVE**

You save lives and we provide an optimal environment. Managing a patient's vital functions is a priority in critical care. LINET ICU beds and mattresses bring advanced solutions to support ICU tasks during CPR, central venous catheter (CVC) insertion, respiratory system management and more.



One-step CPR. Both bed and mattress are ready for cardiopulmonary resuscitation. Quick and safe.



Automatic lateral therapy can help to prevent ARDS, VAP and atelectasis.



In-bed X-ray examination can help to monitor the central venous catheter position.

Multicare

Eleganza 3XC

- ✓ Smart CPR feature
- ✓ X-ray imaging

Eleganza 5

- ✓ Smart CPR feature
- ✓ X-ray imaging

- ✓ Smart CPR feature
- ✓ X-ray imaging
- ✓ Automatic lateral therapy

ICU SAFE

We know how to prevent intensive care-related risks! Our solutions will help you with respiratory complications, falls and pressure ulcers.



OptiCare smart mattress automatically optimises the pressure with respect to individual patient needs.



ErgoFrame® extends the space in the pelvic area and reduces the risk of pressure ulcer development. The position of the cardiac chair is comfortable for the patient while supporting respiration.



Combination of siderails with multi-zone bed exit alarm and i-Brake helps to minimise the risk of falls.

Eleganza 3XC

- ✓ Virtuoso alternating pressure mattress
- ✓ ErgoFrame® positioning feature
- ✓ Cardiac Chair position
- ✓ Siderails with Protector
- ✓ Multi-zone bed exit alarm
- ✓ Brake alarm

Eleganza 5

- ✓ Virtuoso alternating pressure mattress
- ✓ OptiCare smart integrated mattress
- ✓ ErgoFrame® positioning feature
- ✓ Lateral tilt
- ✓ Orthopnoeic Chair position
- ✓ 30°/45° backrest stop
- ✓ Lateral tilt 15°
- ✓ Siderails with Protector
- ✓ Multi-zone bed exit alarm

Multicare

- ✓ Virtuoso alternating pressure mattress
- ✓ OptiCare smart integrated mattress
- ✓ ErgoFrame® positioning feature
- ✓ Automatic lateral therapy
- ✓ Orthopnoeic Chair position
- ✓ 30°/45° backrest stop
- ✓ Automatic lateral therapy
- ✓ Siderails with Protector
- ✓ Multi-zone bed exit alarm
- ✓ i-Brake

icu recovered

Immobilisation syndrome is a barrier to ICU patient recovery. Early mobilisation using innovative ICU bed features safely and effectively accelerates patient recovery.



Mobi-Lift® handle facilitates the sit to stand procedures during mobilisation.



Lateral tilt and siderail override features increase the safety and efficiency of mobilisation.

Mobi-Lift® handle to facilitate getting up from the bed.

Multicare

Eleganza 3XC

- ✓ Mobi-Lift® handle

Eleganza 5

- ✓ Mobi-Lift® handle +
MobiRails +
Lateral tilt for
mobilisation

- ✓ Mobi-Lift® handle +
MobiRails +
Lateral tilt for
mobilisation

icu well

We help to protect ICU staff!
Patient handling and transport
can be performed with
minimum risk
and physical effort.



i-Drive Power – The smart accelerator optimises the drive speed in different drive situations.



The motorised in-bed transport i-Drive Power requires minimum manpower to drive.



The nurse will carry out routine nursing procedures ergonomically and effortlessly.

Eleganza 3XC

- ✓ Scales
- ✓ Foot controls with examination position

Eleganza 5

- ✓ Scales
- ✓ Foot controls with examination position
- ✓ i-Drive Power
- ✓ Lateral tilt

Multicare

- ✓ Scales
- ✓ Foot controls with examination position
- ✓ i-Drive Power
- ✓ Lateral tilt

ICU complete solution

Patient

- Fewer pulmonary complications thanks to ALT®
- Quicker recovery thanks to Mobi-Lift®
- Reduced falls thanks to the siderail concept, Multi-zone bed exit alarm and i-Brake




- Less physical strain and easier patient handling due to the lateral tilt
- Safer and quicker transport thanks to the i-Drive Power
- Easier repositioning due to the Ergoframe®



Staff

Management

- Reduced of extra cost for VAP, falls, infections and pressure ulcers
- Saving staff time for qualified jobs
- Effective management equipment = cost-saving products



 **wissner-
bosserhoff**



Members of LINET Group

LINET spol. s r. o.

Želečnice 5 | 274 01 Slaný | Czech Republic

tel.: +420 312 576 400 | fax: +420 312 522 668 | e-mail: info@linet.com | www.linet.com

